


ALTOS-XL Drive PTA Balloon



- Large Vessel Dilatation
- High Pressure Dilatation
- Low Profile & Short Shoulders

**THE NEW ALTOS-XL DRIVE PTA BALLOON FEATURES
SHORT SHOULDERS & RADIOPAQUE TIP !!!**

Andratec 

Advanced Interventional Technologies



- The ALTOS-XL Drive PTA Balloon is designed for large lumen dilatation. Unique Refolding technology assists the re-insertion.
- The high-pressure performance combined with low profile makes the ALTOS-XL Drive PTA Balloon your choice to treat difficult large vascular and non-vascular lesions.
- Advanced **PET** Material (APM) makes the ALTOS-XL Drive PTA balloon puncture resistant and reliable for multiple dilatations.
- Increased Tip Visibility leading to more accurate Positioning.

Technical Information:

Balloon Material:	NPM New PET Material – 3 to 4 Folds
Balloon Compliance:	Non-compliant
Shaft Diameter:	7F (14&16&18) 8F (20 to 26) 9F (28&30)
Profile:	Refolding Technology
Nominal Pressure:	4 / 4 / 3 / 2 bar
Rated Burst Pressure:	10 / 8 / 6 / 5 / 4 bar
Shaft Length:	125 cm
Recommended Guide Wire	0.035" (0.89 mm)
Recommended Sheath Introducer:	Ranging 7F – 12F

A choice of various sizes (diameters and lengths) will assist your patients needs. Future sizes will be available soon.

ALTOS-XL <i>Drive</i> PTA Catheter Order Codes						
Balloon Ø sheath compatibility (F)	Catheter Shaft	Nom & Rated Bar	20 mm length 125cm	30 mm length 125cm	40 mm length 125cm	60 mm length 125cm
14 mm / 7F	7F	4 / 10 atm	AD-14-20-125	AD-14-30-125	AD-14-40-125	AD-14-60-125
16 mm / 8F	7F	4 / 8 atm	AD-16-20-125	AD-16-30-125	AD-16-40-125	AD-16-60-125
18 mm / 8F	7F	4 / 8 atm	AD-18-20-125	AD-18-30-125	AD-18-40-125	AD-18-60-125
20 mm / 9F	8F	4 / 6 atm	AD-20-20-125	AD-20-30-125	AD-20-40-125	AD-20-60-125
22 mm / 10F	8F	4 / 6 atm	AD-22-20-125	AD-22-30-125	AD-22-40-125	AD-22-60-125
24 mm / 10F	8F	3 / 5 atm	AD-24-20-125	AD-24-30-125	AD-24-40-125	AD-24-60-125
26 mm / 12F	8F	3 / 5 atm	AD-26-20-125	AD-26-30-125	AD-26-40-125	AD-26-60-125
28 mm / 12F	9F	3 / 4 atm	AD-28-20-125		AD-28-40-125	
30 mm / 12F	9F	2 / 4 atm	AD-30-20-125		AD-30-40-125	



Distributed By:



AndraTec GmbH
 Simmernerstr. 70
 D-56075 Koblenz
 Tel: +49 261 29675 172
 Fax: +49 261 29675 173
 www.andratec.com