


Exeter Retrieval Snare




- Unique Snare-Loop Design
- Double Capture Technology
- High Visibility
- Superior Quality

Exeter
Advanced Retrieval Technology

AndraTec *Advanced Interventional Technologies*



1. Atraumatic Nitinol loop with gold-plated tungsten for highest visibility
2. Pre-shaped MP Tip eases navigation
3. Shapable introducer to ease navigation & gain performances
4. **AndraTec** Supportive sleeve avoids kinking when manipulating
5. Luer Connector
6. Extra XXL torque device to ease manipulation
7. One piece kink-resistant wire enhances 1:1 torque for better manipulation
8. Radiopaque platinum inner Marker Ring ensures nearest positioning of introducer



The **Exeter Snare & Exeter Snare Micro** have been especially designed to meet the high requirements during extraction.

The **Exeter Retrieval Technology** is your first choice safety device when it comes to successful retrieval.

Exeter
Advanced Retrieval Technology

Exeter Snare Order Codes

Product Code	Diameter (mm)	Length (cm)	Introducer Sheath (cm)	Introducer Sheath (F)
ES - 5	5mm (3-5mm vessel)	125	110	4
ES-10	10mm (6-10mm vessel)	125	110	4
ES-15	15mm (11-15mm vessel)	125	110	4
ES-20	20mm (16-20mm vessel)	125	110	4
ES-25	25mm (21-25mm vessel)	125	110	5
ES-30	30mm (26-30mm vessel)	125	110	5
ES-35	35mm (31-35mm vessel)	125	110	5
Product Code	Diameter (mm)	Length (cm)	Introducer Sheath (cm)	Introducer Sheath (F)
ES - 25/65	25mm (21-25mm vessel)	65	50	5
ES-25/150	25mm (21-25mm vessel)	150	135	5
Product Code	Diameter (mm)	Length (cm)	Introducer Sheath (cm)	Introducer Sheath (F)
ESM-2	2mm (1-2mm vessel)	175	150	3-2.3
ESM-4	4mm (3-4mm vessel)	175	150	3-2.3
ESM-7	7mm (5-7mm vessel)	175	150	3-2.3

CE 0481

Distributed By:

AndraTec
Advanced Interventional Technologies

AndraTec GmbH
Simmernerstr. 70
D-56075 Koblenz
Tel: +49 261 29675 172
Fax: +49 261 29675 173
www.andratec.com